

## ProteoAnnotator user guide

- 1- Download the latest version from: <http://www.proteoannotator.org/?q=installation>
- 2- Extract the zip file. Example, extract ProteoAnnotator-1.0.zip
- 3- Download the test files from: <http://www.proteoannotator.org/?q=installation>
- 4- Extract the zip file. Example, extract the ProteoAnnotator\_test\_files.zip
- 5- Change the directory to the extracted ProteoAnnotator-1.0.zip.
- 6- Run the following command:

```
java -Xmx10G -jar mzidentml-lib-1.6.10-SNAPSHOT.jar ProteoAnnotator -inputGFF
PATH_TO_ProteoAnnotator_test_files\canonical_gene_model\ToxoDB-10.0_TgondiiME49.gff -
spectrum_files PATH_TO_ProteoAnnotator_test_files\mgf\Orbitap-1DE\combined.mgf -
outputFolder PATH_TO_ProteoAnnotator_test_files\output_results -inputPredicted
"PATH_TO_ProteoAnnotator_test_files\non_canonical_gene_model\augustus\augustus.gff;PAT
H_TO_ProteoAnnotator_test_files\non_canonical_gene_model\augustus\augustus.aa##PATH_TO
_ProteoAnnotator_test_files\non_canonical_gene_model\glimmer\Toxo-Glimmer-ME49.gff" -
prefix orbitap-1DE toxo10 -compress false -searchParameters
"PATH_TO_ProteoAnnotator_test_files\search.txt" -peptideThreshValue 0.01 -
proteinThreshValue 0.01 -enableMsgf true
```

**IMPORTANT:** Replace **PATH\_TO\_ProteoAnnotator\_test\_files** with the actual path to the unzipped ProteoAnnotator\_test\_files folder.

### Example (change it to match your folder structures):

```
java -Xmx10G -jar mzidentml-lib-1.6.10-SNAPSHOT.jar ProteoAnnotator -inputGFF
C:\ProteoAnnotator_test_files\canonical_gene_model\ToxoDB-10.0_TgondiiME49.gff -
spectrum_files C:\ProteoAnnotator_test_files\mgf\Orbitap-1DE\combined.mgf -
outputFolder C:\ProteoAnnotator_test_files\output_results -inputPredicted
"C:\ProteoAnnotator_test_files\non_canonical_gene_model\augustus\augustus.gff;C:\Prote
oAnnotator_test_files\non_canonical_gene_model\augustus\augustus.aa##C:\ProteoAnnotato
r_test_files\non_canonical_gene_model\glimmer\Toxo-Glimmer-ME49.gff" -prefix orbitap-
1DE_toxo10_ -compress false -searchParameters
"C:\ProteoAnnotator_test_files\search.txt" -peptideThreshValue 0.01 -
proteinThreshValue 0.01 -enableMsgf true
```

In this example the "**PATH\_TO\_ProteoAnnotator\_test\_files**" is "C:\ProteoAnnotator\_test\_files"

Parameter	optional/ mandatory	Explanation
<b>-prefix</b>	Optional	A prefix to be attached to the output file names
<b>-inputGFF_A</b>	Mandatory	The canonical GFF file
<b>-inputFasta_A</b>	Optional if the canonical GFF contains the FASTA	The protein database
<b>-outputFolder</b>	Mandatory	The output folder for the analysis
<b>-spectrum_files</b>	Mandatory	The MGF file ( <b>only one input file currently supported, later builds will incorporate multiple files to be merged</b> ).
<b>-searchParameters</b>	Mandatory	The search parameters file to be used for the search.
<b>-inputPredicted</b>	Optional	The non-canonical gene models, these are a set of GFF/FASTA files. The GFF and FASTA are separated by ';' and the pairs are separated by

		'# #'
- peptideThreshValue	Mandatory	Peptide FDR Threshold value
- proteinThreshValue	Mandatory	Protein FDR Threshold value
-enableMsgf	Mandatory	Enable/Disable MSGF+